

# **ConvertAble<sup>™</sup> Series**

A brand new concept in nasal masks



More Convenience, Choice, Cost Savings

# **ConvertAble<sup>™</sup> Series**

## The choice of two mask cushions

The ConvertAble<sup>™</sup> Series is a ResMed exclusive, offering the choice of two highperformance cushions—Mirage<sup>™</sup> SoftGel or Mirage Activa<sup>™</sup> LT—on one durable frame. The benefits:

- Convenience You can easily swap mask cushions without costly replacements
- Choice Two mask options empower clinicians to choose the best seal for each patient
- **Cost savings** Common parts streamline inventory and interchangeable cushions lower switching costs

There's less inventory, less hassle and more opportunity to drive replacement. Only ResMed supplies the ConvertAble frame—providing the best possible solutions for your patients.

### All backed by ResMed's guarantee of quality and durability.

### Mirage Activa<sup>™</sup> LT

NASAL MASK

### Activa LT with ActiveCell<sup>™</sup> for active sleepers

- ActiveCell inflatable chamber provides an active seal to move with the patient as they sleep
- Functions as an independent suspension helping to keep the cushion in place
- Thick inner cushion provides comfortable support and stability
- Thin outer cushion seals lightly against the face and prevents facial marking





### on one durable frame.

#### ResMed's exclusive ConvertAble frame

MED

- MicroFit dial 24 positions for easy and precise adjustment
- Flexible forehead pads provide support and comfort
- One-snap elbow 360 degree rotation for flexibility
- Headgear clips quick-release for added convenience
- Diffused vents direct airflow quietly away for comfort
- Quiet provide peaceful sleep for patient and bed partner

### Mirage<sup>™</sup> SoftGel

### NASAL MASK

#### SoftGel comfort for strong compliance

- · Easy, intuitive first-time fit
- DoubleGel<sup>™</sup> cushion automatically conforms to patient's face, offering a personalized fit and seal
- Two layers for comfort and stability
- Patient-preferred over competitor masks

**Softness:** The clear layer lightly seals and conforms to unique facial features

**Support:** The blue gel layer provides support and stability





### **PRODUCT CODES**

### ConvertAble Pack Mirage SoftGel & Mirage Activa LT Cushion

Small	61604
Medium	61609
Large	61615
Large-Wide	61620

### Mirage SoftGel Nasal Mask

Small	61600
Medium	61601
Large	61602
Large-Wide	61603

### Mirage Activa LT Nasal Mask

Small	60182
Medium	60148
Large	60149
Large-Wide	60150

#### **REIMBURSEMENT CODES**

Nasal mask use with positive airway pressure device, 1 per 3 months	A7034
Replacement cushion for nasal mask, 2 per month	A7032
Headgear used with positive airway pressure device, 1 per 6 months	A7035

Example of a CMS replacement schedule that meets current Medicare quidelines, which allows replacement of up to two nasal cushions per month.

Month 1	Month 2	Month 3
A7034 A7032 A7135		A7032
Month 4	Month 5	Month 6
A7032	A7032	A7032
Month 7	Month 8	Month 9
A7034 A7032 A7135		A7032
Month 10	Month 11	Month 12
A7032	A7032	A7032

The Medicare schedule does not state how often Medicare recipients should replace accessories, but it does state how often they can be replaced if necessary. It is always the provider's responsibility to determine and submit appropriate codes, charges, modifiers and bills for the services that were rendered. Before filing any claims, providers should verify current requirements and policies with the payor.

\*Either SoftGel or Activa LT cushion can be dispensed on an ongoing basis.

ResMed Corp San Diego, CA, USA +1 858 836 5000 or 1 800 424 0737 (toll free). ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). See www.resmed.com for other ResMed locations worldwide. ConvertAble, DoubleGel, and Mirage Activa are trademarks of ResMed Ltd and registered in the U.S. Patent and Trademark Office. ©2009 ResMed Ltd. Specifications may change without notice. 1013511/3 09 11