

Our Mojo® Vented and Non-Vented masks feature a patented AIR°gel cushion, flexible shell with Custom Fit Technology™ and an Adjustable Tension Wheel. Together, they provide an optimal fit with maximum comfort. Fine-tuned fit and function — have it all in a full-face mask.

Your Mask, Your Way

AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.

Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

Flex, fit, forget.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

AIR°gel Forehead Cushion


Adjustable Tension Wheel

Removes pressure from the bridge of the nose for a custom, comfortable fit.

Easy-Release Magnetic Headgear Clips

Lightweight, breathable headgear disconnects from mask easily with one hand.

Vented and Non-Vented Versions Available



Made in the USA

What makes our Mojo Mask different?

- Custom Fit Technology is unique to Sleepnet®. Our soft-shell mask has a shapeable wire allowing for optimal flexibility, fit and comfort.
- The Adjustable Tension Wheel is exclusive to the Sleepnet Mojo. A wide, continuous range of adjustment ensures a custom fit and removes pressure from the bridge of the nose.
- Available in three sizes with magnetic connectors and simple, breathable quick-release headgear.

About AIR°gel

Patented AIR°gel cushion delivers an effective seal while minimizing pressure points on the skin.



How does Mojo stand apart from the competition?

Features	Sleepnet Mojo Mask	ResMed Mirage Quattro™ Fx	Fisher & Paykel Forma Mask	Philips Respironics Comfort Gel Full
Custom Fit Technology	●			
Patented AIR°gel	●			
Fine tune forehead adjustment	●			
Easy-release magnetic headgear clips	●			
Available in 3 sizes	●	●	●	●
100% made in America	●			

Mojo Masks with Headgear Specifications:

Product Number	Mask Style	Mask Size	Mask Weight	Mask Dimensions
50835	Vented	Small	144 g	3.73"D x 3.78"W x 6.50"H
50836	Vented	Medium	148 g	3.72"D x 3.76"W x 6.75"H
50837	Vented	Large	156 g	3.79"D x 3.76"W x 7.00"H
50920	Non-Vented	Small	144 g	3.73"D x 3.78"W x 6.50"H
50921	Non-Vented	Medium	148 g	3.72"D x 3.76"W x 6.75"H
50922	Non-Vented	Large	156 g	3.79"D x 3.76"W x 7.00"H

CE 0123



Does not contain natural rubber latex

iQ® Blue™

The Smart Choice for an Off-the-Shelf Custom Fit Nasal Mask.



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646
Visit www.sleepnetmasks.com to learn more.
Sleepnet Corporation, 5 Merrill Industrial Drive,
Hampton, New Hampshire 03842 USA

